

SNX1 Rabbit mAb

Catalog No: #59870

Package Size: #59870-1 50ul #59870-2 100ul


Description

Product Name	SNX1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	SNX1 Antibody detects endogenous levels of total SNX1
Immunogen Description	A synthesized peptide derived from human SNX1
Other Names	Snx1; SNX1A; Sorting nexin 1; Sorting nexin 1A; Vps5;
Accession No.	Uniprot:Q13596
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of SNX1 expression in 293T cell lysate.

Background

May be involved in several stages of intracellular trafficking. Plays a role in targeting ligand-activated EGFR to the lysosomes for degradation after endocytosis from the cell surface and release from the Golgi.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.