

VAMP1 Rabbit mAb

Catalog No: #59873

Package Size: #59873-1 50ul #59873-2 100ul


Description

Product Name	VAMP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	VAMP1 Antibody detects endogenous levels of total VAMP1
Immunogen Description	A synthesized peptide derived from human VAMP1
Other Names	SYB1; Synaptobrevin 1; Vamp1;
Accession No.	Uniprot:P23763
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of VAMP1 expression in Human cerebellum lysate.

Background

Involved in the targeting and/or fusion of transport vesicles to their target membrane.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.