

NUP153 Rabbit mAb

Catalog No: #59879

Package Size: #59879-1 50ul #59879-2 100ul


Description

Product Name	NUP153 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	Nup153 Antibody detects endogenous levels of total Nup153
Immunogen Description	A synthesized peptide derived from human Nup153
Other Names	HNUP153; Nup153;
Accession No.	Uniprot:P49790
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Nup153 expression in K562 cell lysate.

Background

Possible DNA-binding subunit of the nuclear pore complex (NPC). The repeat-containing domain may be involved in anchoring components of the pore complex to the pore membrane.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.