

E2F4 Rabbit mAb

Catalog No: #59891

Package Size: #59891-1 50ul #59891-2 100ul


Description

Product Name	E2F4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	E2F4 Antibody detects endogenous levels of total E2F4
Immunogen Description	A synthesized peptide derived from human E2F4
Other Names	E2F 4; E2F transcription factor 4; E2F transcription factor 4 p107/p130 binding; E2F4;
Accession No.	Uniprot:Q16254
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of E2F4 expression in Raji cell lysate.

Background

Transcription activator that binds DNA cooperatively with DP proteins through the E2 recognition site, 5'-TTTC[CG]CGC-3' found in the promoter region of a number of genes whose products are involved in cell cycle regulation or in DNA replication.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.