

GST3 / GSTP1 Rabbit mAb

Catalog No: #59892

Package Size: #59892-1 50ul #59892-2 100ul


Description

Product Name	GST3 / GSTP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	GST3 / GST pi Antibody detects endogenous levels of total GST3 / GST pi
Immunogen Description	A synthesized peptide derived from human GST3 / GST pi
Other Names	Deafness; DFN7; FAEES3; GST3; GSTP; Gstp1; PI; X linked 7;
Accession No.	Uniprot:P09211
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of GST3 / GST pi expression in PC3 cell lysate.

Background

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.