

## CD27 Rabbit mAb

Catalog No: #59901

Package Size: #59901-1 50ul #59901-2 100ul


## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | CD27 Rabbit mAb                                                                                    |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB ICC/IF FC                                                                                       |
| Species Reactivity    | Human Rat                                                                                          |
| Specificity           | CD27 Antibody detects endogenous levels of total CD27                                              |
| Immunogen Description | A synthesized peptide derived from human CD27                                                      |
| Other Names           | CD27; CD27L receptor; LPFS2; S152; T14; TNFRSF7; TNFSF7; Tp55;                                     |
| Accession No.         | Uniprot:P26842                                                                                     |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                     |

## Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

## Images



Western blot analysis of CD27 expression in Ramos cell lysate.

## Background

Receptor for CD70/CD27L. May play a role in survival of activated T-cells. May play a role in apoptosis through association with SIVA1.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.