

CXCR5 Rabbit mAb

Catalog No: #59904

Package Size: #59904-1 50ul #59904-2 100ul


Description

Product Name	CXCR5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	CXCR5 Antibody detects endogenous levels of total CXCR5
Immunogen Description	A synthesized peptide derived from human CXCR5
Other Names	BLR1; CD185; Cxcr5; MDR15;
Accession No.	Uniprot:P32302
Concentration	0.4mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:1000 - 1:2000; IF/ICC 1:50 - 1:200; FC 1:20-1:100

Images



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Background

Cytokine receptor that binds to B-lymphocyte chemoattractant (BLC). Involved in B-cell migration into B-cell follicles of spleen and Peyer patches but not into those of mesenteric or peripheral lymph nodes. May have a regulatory function in Burkitt lymphoma (BL) lymphomagenesis and/or B-cell differentiation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.