

BAF57/SMARCE1 Rabbit mAb

Catalog No: #59905

Package Size: #59905-1 50ul #59905-2 100ul


Description

Product Name	BAF57/SMARCE1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	BAF57/SMARCE1 Antibody detects endogenous levels of total BAF57/SMARCE1
Immunogen Description	A synthesized peptide derived from human BAF57/SMARCE1
Other Names	BAF57; SMARCE1;
Accession No.	Uniprot:Q969G3
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of BAF57/SMARCE1 expression in Raji cell lysate.

Background

Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-nucleosome topology). Belongs to the neural progenitors-specific chromatin remodeling complex (npBAF complex) and the neuron-specific chromatin remodeling complex (nBAF complex).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.