

PEN2 Rabbit mAb

Catalog No: #59910

Package Size: #59910-1 50ul #59910-2 100ul


Description

Product Name	PEN2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	PEN2 Antibody detects endogenous levels of total PEN2
Immunogen Description	A synthesized peptide derived from human PEN2
Other Names	PEN 2; PSEN2; psenen;
Accession No.	Uniprot:Q9NZ42
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of PEN2 expression in Human placenta lysate.

Background

Essential subunit of the gamma-secretase complex, an endoprotease complex that catalyzes the intramembrane cleavage of integral membrane proteins such as Notch receptors and APP (beta-amyloid precursor protein).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.