

Mad2L1 Rabbit mAb

Catalog No: #59915

Package Size: #59915-1 50ul #59915-2 100ul


Description

Product Name	Mad2L1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Mad2L1 Antibody detects endogenous levels of total Mad2L1
Immunogen Description	A synthesized peptide derived from human Mad2L1
Other Names	HsMAD2; MAD 2; MAD 2A; MAD2 like 1; Mad2l1;
Accession No.	Uniprot:Q13257
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of Mad2L1 expression in A431 cell lysate.

Background

Component of the spindle-assembly checkpoint that prevents the onset of anaphase until all chromosomes are properly aligned at the metaphase plate. Required for the execution of the mitotic checkpoint which monitors the process of kinetochore-spindle attachment and inhibits the activity of the anaphase promoting complex by sequestering CDC20 until all chromosomes are aligned at the metaphase plate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.