

MBD3 Rabbit mAb

Catalog No: #59917

Package Size: #59917-1 50ul #59917-2 100ul


Description

Product Name	MBD3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	MBD3 Antibody detects endogenous levels of total MBD3
Immunogen Description	A synthesized peptide derived from human MBD3
Other Names	Mbd3;
Accession No.	Uniprot:O95983
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

Images



Western blot analysis of MBD3 expression in 293T cell lysate.

Background

Does not bind DNA by itself. Recruits histone deacetylases and DNA methyltransferases. Acts as transcriptional repressor and plays a role in gene silencing.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.