

## FYB Rabbit mAb

Catalog No: #59937

Package Size: #59937-1 50ul #59937-2 100ul


## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | FYB Rabbit mAb                                                                                     |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB IHC ICC/IF                                                                                      |
| Species Reactivity    | Mouse, Human                                                                                       |
| Specificity           | FYB Antibody detects endogenous levels of total FYB                                                |
| Immunogen Description | A synthesized peptide derived from human FYB                                                       |
| Other Names           | ADAP; Fyb; FYN binding protein; p120/p130; SLAP130;                                                |
| Accession No.         | Uniprot:O15117                                                                                     |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                     |

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

## Images



Western blot analysis of FYB expression in TF-1 cell lysate.

## Background

Acts as an adapter protein of the FYN and LCP2 signaling cascades in T cells. Modulates the expression of interleukin-2 (IL-2). Involved in platelet activation. Prevents the degradation of SKAP1 and SKAP2.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.