

IFNGR1 Rabbit mAb

Catalog No: #59954

Package Size: #59954-1 50ul #59954-2 100ul


Description

Product Name	IFNGR1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	IFNGR1 Antibody detects endogenous levels of total IFNGR1
Immunogen Description	A synthesized peptide derived from human IFNGR1
Other Names	CD119; CDw119; IFN gamma R alpha; IFNGR1;
Accession No.	Uniprot:P15260
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:500

Images



Western blot analysis of IFNGR1 expression in MCF7 cell lysate.

Background

Receptor for interferon gamma. Two receptors bind one interferon gamma dimer.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.