

ACADM Rabbit mAb

Catalog No: #59961

Package Size: #59961-1 50ul #59961-2 100ul


Description

Product Name	ACADM Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC IP
Species Reactivity	Human Mouse Rat
Specificity	ACADM Antibody detects endogenous levels of total ACADM
Immunogen Description	A synthesized peptide derived from human ACADM
Other Names	ACAD1; Acadm; MCAD; MCADH;
Accession No.	Uniprot:P11310
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of ACADM expression in HepG2 cell lysate.

Background

This enzyme is specific for acyl chain lengths of 4 to 16.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.