

TEAD1 Rabbit mAb

Catalog No: #60024

Package Size: #60024-1 50ul #60024-2 100ul


Description

Product Name	TEAD1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	TEAD1 Antibody detects endogenous levels of total TEAD1
Immunogen Description	A synthesized peptide derived from human TEAD1
Other Names	NTEF1; REF1; TCF13; TEA domain family member 1; TEAD1; TEF 1; Transcription factor 13;
Accession No.	Uniprot:P28347
Calculated MW	52kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200

Images



Western blot analysis of TEAD1 expression in 293T cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Background

Transcription factor which plays a key role in the Hippo signaling pathway, a pathway involved in organ size control and tumor suppression by restricting proliferation and promoting apoptosis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.