

Inhibin beta A Rabbit mAb

Catalog No: #60025

Package Size: #60025-1 50ul #60025-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Inhibin beta A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Inhibin beta A Antibody detects endogenous levels of total Inhibin beta A
Immunogen Description	A synthesized peptide derived from human Inhibin beta A
Other Names	EDF; FRP; INHBA;
Accession No.	Uniprot:P08476
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of Inhibin beta A expression in Human fetal brain lysate.

Background

Inhibins and activins inhibit and activate, respectively, the secretion of follitropin by the pituitary gland. Inhibins/activins are involved in regulating a number of diverse functions such as hypothalamic and pituitary hormone secretion, gonadal hormone secretion, germ cell development and maturation, erythroid differentiation, insulin secretion, nerve cell survival, embryonic axial development or bone growth, depending on their subunit composition.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.