

PLA2G2A Rabbit mAb

Catalog No: #60040

Package Size: #60040-1 50ul #60040-2 100ul


Description

Product Name	PLA2G2A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	PLA2G2A Antibody detects endogenous levels of total PLA2G2A
Immunogen Description	A synthesized peptide derived from human PLA2G2A
Other Names	MOM1; PLA2; PLA2B; PLA2L; PLA2S; PLAS1; sPLA2;
Accession No.	Uniprot:P14555
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of PLA2G2A expression in A549 cell lysate.

Background

Thought to participate in the regulation of the phospholipid metabolism in biomembranes including eicosanoid biosynthesis. Catalyzes the calcium-dependent hydrolysis of the 2-acyl groups in 3-sn-phosphoglycerides.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.