

Nidogen 1 Rabbit mAb

Catalog No: #60041

Package Size: #60041-1 50ul #60041-2 100ul


Description

Product Name	Nidogen 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Nidogen 1 Antibody detects endogenous levels of total Nidogen 1
Immunogen Description	A synthesized peptide derived from human Nidogen 1
Other Names	Enactin; Ent; Entactin 1; Entactin; NID 1; NID; Nid1; Nidogen 1; Nidogen;
Accession No.	Uniprot:P14543
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of Nidogen 1 expression in Human fetal heart lysate.

Background

Sulfated glycoprotein widely distributed in basement membranes and tightly associated with laminin. Also binds to collagen IV and perlecan. It probably has a role in cell-extracellular matrix interactions.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.