

NME2 Rabbit mAb

Catalog No: #60046

Package Size: #60046-1 50ul #60046-2 100ul


Description

Product Name	NME2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	NME2 Antibody detects endogenous levels of total NME2
Immunogen Description	A synthesized peptide derived from human NME2
Other Names	C myc purine binding transcription factor PUF; MGC2212; NDKB; NDP kinase B; NDPKB; NM23B; Nucleoside diphosphate kinase B; PUF;
Accession No.	Uniprot:P22392
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of NME2 expression in LnCaP cell lysate.

Background

Major role in the synthesis of nucleoside triphosphates other than ATP. Negatively regulates Rho activity by interacting with AKAP13/LBC.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.