

MLKL Rabbit mAb

Catalog No: #60050

Package Size: #60050-1 50ul #60050-2 100ul


Description

Product Name	MLKL Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	MLKL Antibody detects endogenous levels of total MLKL
Immunogen Description	A synthesized peptide derived from MLKL
Other Names	hMLKL; Mixed lineage kinase domain like; Mlkl;
Accession No.	Uniprot:Q8NB16
Calculated MW	54kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000 IP 1:50

Images



Western blot analysis of MLKL expression in Hela cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Background

Pseudokinase that plays a key role in TNF-induced necroptosis, a programmed cell death process.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.