

UBE4B Rabbit mAb

Catalog No: #60061

Package Size: #60061-1 50ul #60061-2 100ul


Description

Product Name	UBE4B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	UBE4B Antibody detects endogenous levels of total UBE4B
Immunogen Description	A synthesized peptide derived from UBE4B
Other Names	HDNB1; Ube4b; UBOX3; UFD2; UFD2a; Ufd2p;
Accession No.	Uniprot:O95155
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of UBE4B expression in (1) 293T cell lysate; (2) NIH/3T3 cell lysate.

Background

Binds to the ubiquitin moieties of preformed conjugates and catalyzes ubiquitin chain assembly in conjunction with E1, E2, and E3.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.