

Desmocollin?3 Rabbit mAb

Catalog No: #60063

Package Size: #60063-1 50ul #60063-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Desmocollin?3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	Desmocollin?3 Antibody detects endogenous levels of total Desmocollin?3
Immunogen Description	A synthesized peptide derived from Desmocollin?3
Other Names	CDHF3; Desmocollin3; Desmocollin4; DSC; DSC1; DSC2; DSC3; DSC4; HTCP;
Accession No.	Uniprot:Q14574
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Desmocollin?3 expression in A431 cell lysate.

Background

Component of intercellular desmosome junctions. Involved in the interaction of plaque proteins and intermediate filaments mediating cell-cell adhesion. May contribute to epidermal cell positioning (stratification) by mediating differential adhesiveness between cells that express different isoforms.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.