

IKZF3 Rabbit mAb

Catalog No: #60087

Package Size: #60087-1 50ul #60087-2 100ul


Description

Product Name	IKZF3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	IKZF3 Antibody detects endogenous levels of total IKZF3
Immunogen Description	A synthesized peptide derived from human IKZF3
Other Names	AIO; IKZF3; ZNFN1A3;
Accession No.	Uniprot:Q9UKT9
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of IKZF3 expression in Raji cell lysate.

Background

Transcription factor that plays an important role in the regulation of lymphocyte differentiation. Plays an essential role in regulation of B-cell differentiation, proliferation and maturation to an effector state. Involved in regulating BCL2 expression and controlling apoptosis in T-cells in an IL2-dependent manner.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.