

Cyclin H/p34 Rabbit mAb

Catalog No: #60106

Package Size: #60106-1 50ul #60106-2 100ul


Description

Product Name	Cyclin H/p34 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Cyclin H/p34 Antibody detects endogenous levels of total Cyclin H/p34
Immunogen Description	A synthesized peptide derived from human Cyclin H/p34
Other Names	CAK; CAK complex subunit; CDK activating kinase; CDK activating kinase complex subunit; CyclinH; MO15 associated protein; p34; p36; p37;
Accession No.	Uniprot:P51946
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of Cyclin H/p34 expression in A431 cell lysate.

Background

Regulates CDK7, the catalytic subunit of the CDK-activating kinase (CAK) enzymatic complex. CAK activates the cyclin-associated kinases CDK1, CDK2, CDK4 and CDK6 by threonine phosphorylation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.