

HMG4 Rabbit mAb

Catalog No: #60114

Package Size: #60114-1 50ul #60114-2 100ul


Description

Product Name	HMG4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse
Specificity	HMG4 Antibody detects endogenous levels of total HMG4
Immunogen Description	A synthesized peptide derived from human HMG4
Other Names	HMG 4; HMG2A; Hmgb3;
Accession No.	Uniprot:O15347
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of HMG4 expression in HeLa cell lysate.

Background

Binds preferentially single-stranded DNA and unwinds double stranded DNA.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.