

Cip4 Rabbit mAb

Catalog No: #60117

Package Size: #60117-1 50ul #60117-2 100ul

Description

Product Name	Cip4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse
Specificity	Cip4 Antibody detects endogenous levels of total Cip4
Immunogen Description	A synthesized peptide derived from human Cip4
Other Names	Cdc42 interacting protein 4; CIP4; DCIP4; hSTP; STOT; STP; Thyroid receptor interacting protein 10; trip10;
Accession No.	Uniprot:Q15642
Calculated MW	75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Background

Required for translocation of GLUT4 to the plasma membrane in response to insulin signaling (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.