

Nodal Rabbit mAb

Catalog No: #60120

Package Size: #60120-1 50ul #60120-2 100ul


Description

Product Name	Nodal Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Nodal Antibody detects endogenous levels of total Nodal
Immunogen Description	A synthesized peptide derived from human Nodal
Other Names	HTX5; Nodal;
Accession No.	Uniprot:Q96S42
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of Nodal expression in SH-SY5Y cell lysate.

Background

Essential for mesoderm formation and axial patterning during embryonic development.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.