

HDLBP Rabbit mAb

Catalog No: #60129

Package Size: #60129-1 50ul #60129-2 100ul


Description

Product Name	HDLBP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	HDLBP Antibody detects endogenous levels of total HDLBP
Immunogen Description	A synthesized peptide derived from human HDLBP
Other Names	HBP; HDL binding protein; HDLBP; VGL;
Accession No.	Uniprot:Q00341
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of HDLBP expression in Jurkat cell lysate.

Background

Appears to play a role in cell sterol metabolism. It may function to protect cells from over-accumulation of cholesterol.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.