

GNA13 Rabbit mAb

Catalog No: #60132

Package Size: #60132-1 50ul #60132-2 100ul


Description

Product Name	GNA13 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	GNA13 Antibody detects endogenous levels of total GNA13
Immunogen Description	A synthesized peptide derived from human GNA13
Other Names	G alpha 13; GNA13;
Accession No.	Uniprot:Q14344
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of GNA13 expression in HepG2 cell lysate.

Background

Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.