

LTK Rabbit mAb

Catalog No: #60135

Package Size: #60135-1 50ul #60135-2 100ul


Description

Product Name	LTK Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	LTK Antibody detects endogenous levels of total LTK
Immunogen Description	A synthesized peptide derived from human LTK
Other Names	Ltk; TYK1;
Accession No.	Uniprot:P29376
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of LTK expression in Raji cell lysate.

Background

Orphan receptor with a tyrosine-protein kinase activity. The exact function of this protein is not known. Studies with chimeric proteins (replacing its extracellular region with that of several known growth factor receptors, such as EGFR and CSFIR) demonstrate its ability to promote growth and specifically neurite outgrowth, and cell survival.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.