

RNF14 Rabbit mAb

Catalog No: #60139

Package Size: #60139-1 50ul #60139-2 100ul


Description

Product Name	RNF14 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	RNF14 Antibody detects endogenous levels of total RNF14
Immunogen Description	A synthesized peptide derived from human RNF14
Other Names	ARA54; HFB30; Rnf14; TRIAD2;
Accession No.	Uniprot:Q9UBS8
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of RNF14 expression in Jurkat cell lysate.

Background

Might act as an E3 ubiquitin-protein ligase which accepts ubiquitin from specific E2 ubiquitin-conjugating enzymes and then transfers it to substrates, which could be nuclear proteins. Could play a role as a coactivator for androgen- and, to a lesser extent, progesterone-dependent transcription.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.