

PLSCR3 Rabbit mAb

Catalog No: #60141

Package Size: #60141-1 50ul #60141-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PLSCR3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	PLSCR3 Antibody detects endogenous levels of total PLSCR3
Immunogen Description	A synthesized peptide derived from human PLSCR3
Other Names	PLS3; Plscr3;
Accession No.	Uniprot:Q9NRY6
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of PLSCR3 expression in BxPC 3 cell lysate.

Background

May mediate accelerated ATP-independent bidirectional transbilayer migration of phospholipids upon binding calcium ions that results in a loss of phospholipid asymmetry in the plasma membrane. May play a central role in the initiation of fibrin clot formation, in the activation of mast cells and in the recognition of apoptotic and injured cells by the reticuloendothelial system.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.