

DNAJC12 Rabbit mAb

Catalog No: #60145

Package Size: #60145-1 50ul #60145-2 100ul

Description

Product Name	DNAJC12 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse
Specificity	DNAJC12 Antibody detects endogenous levels of total DNAJC12
Immunogen Description	A synthesized peptide derived from human DNAJC12
Other Names	DNAJC12; HPANBH4; JDP1; mJDP1;
Accession No.	Uniprot:Q9UKB3
Calculated MW	23kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Background

Contains 1 J domain.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.