

KIF5B Rabbit mAb

Catalog No: #60149

Package Size: #60149-1 50ul #60149-2 100ul

Description

Product Name	KIF5B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	KIF5B Antibody detects endogenous levels of total KIF5B
Immunogen Description	A synthesized peptide derived from human KIF5B
Other Names	KIF5B; Kinesin1; KINH; KNS; KNS1; UKHC;
Accession No.	Uniprot:P33176
Calculated MW	120kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Background

Microtubule-dependent motor required for normal distribution of mitochondria and lysosomes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.