

## MRPS31 Rabbit mAb

Catalog No: #60153

Package Size: #60153-1 50ul #60153-2 100ul


## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | MRPS31 Rabbit mAb                                                                                  |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB                                                                                                 |
| Species Reactivity    | Human Mouse Rat                                                                                    |
| Specificity           | MRPS31 Antibody detects endogenous levels of total MRPS31                                          |
| Immunogen Description | A synthesized peptide derived from human MRPS31                                                    |
| Other Names           | IMOGN38; MRPS31; S31mt;                                                                            |
| Accession No.         | Uniprot:Q92665                                                                                     |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                     |

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of MRPS31 expression in HepG2 cell lysate.

## Background

Component of the mitochondrial ribosome small subunit (28S) which comprises a 12S rRNA and about 30 distinct proteins.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.