

C3orf38 Rabbit mAb

Catalog No: #60155

Package Size: #60155-1 50ul #60155-2 100ul

Description

Product Name	C3orf38 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	C3orf38 Antibody detects endogenous levels of total C3orf38
Immunogen Description	A synthesized peptide derived from human C3orf38
Other Names	C3orf38;
Accession No.	Uniprot:Q5JPI3
Calculated MW	37kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Background

May be involved in apoptosis regulation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.