

ACAT1 Rabbit mAb

Catalog No: #60162

Package Size: #60162-1 50ul #60162-2 100ul

Description

Product Name	ACAT1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	ACAT1 Antibody detects endogenous levels of total ACAT1
Immunogen Description	A synthesized peptide derived from human ACAT1
Other Names	ACAT; acat1; MAT; RATACAL; THIL;
Accession No.	Uniprot:P24752
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

Background

Plays a major role in ketone body metabolism.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.