

PCTAIRE1 Rabbit mAb

Catalog No: #60171

Package Size: #60171-1 50ul #60171-2 100ul

Description

Product Name	PCTAIRE1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	PCTAIRE1 Antibody detects endogenous levels of total PCTAIRE1
Immunogen Description	A synthesized peptide derived from human PCTAIRE1
Other Names	Cdk16; PCTAIRE; PCTAIRE1; PCTGAIRE; PCTK1;
Accession No.	Uniprot:Q00536
Calculated MW	62-68kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Background

May play a role in signal transduction cascades in terminally differentiated cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.