

AAMP Rabbit mAb

Catalog No: #60181

Package Size: #60181-1 50ul #60181-2 100ul

Description

Product Name	AAMP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	AAMP Antibody detects endogenous levels of total AAMP
Immunogen Description	A synthesized peptide derived from human AAMP
Other Names	AAMP;
Accession No.	Uniprot:Q13685
Calculated MW	52kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Background

Plays a role in angiogenesis and cell migration. In smooth muscle cell migration, may act through the RhoA pathway.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.