

UFD1L Rabbit mAb

Catalog No: #60188

Package Size: #60188-1 50ul #60188-2 100ul


Description

Product Name	UFD1L Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	UFD1L Antibody detects endogenous levels of total UFD1L
Immunogen Description	A synthesized peptide derived from human UFD1L
Other Names	UB fusion protein 1; UFD1; UFD1L;
Accession No.	Uniprot:Q92890
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of UFD1L expression in HeLa cell lysate.

Background

Essential component of the ubiquitin-dependent proteolytic pathway which degrades ubiquitin fusion proteins. The ternary complex containing UFD1L, VCP and NPLOC4 binds ubiquitinated proteins and is necessary for the export of misfolded proteins from the ER to the cytoplasm, where they are degraded by the proteasome.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.