

PPM1G Rabbit mAb

Catalog No: #60196

Package Size: #60196-1 50ul #60196-2 100ul


Description

Product Name	PPM1G Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Rat
Specificity	PPM1G Antibody detects endogenous levels of total PPM1G
Immunogen Description	A synthesized peptide derived from human PPM1G
Other Names	PP2C; PP2C gamma; PP2CG; PPM1G; PPP2CG; Protein phosphatase 1C; Protein phosphatase 1G;
Accession No.	Uniprot:O15355
Calculated MW	75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of PPM1G expression in Jurkat cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Background

Belongs to the PP2C family.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.