

IQGAP2 Rabbit mAb

Catalog No: #60197

Package Size: #60197-1 50ul #60197-2 100ul

Description

Product Name	IQGAP2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	IQGAP2 Antibody detects endogenous levels of total IQGAP2
Immunogen Description	A synthesized peptide derived from human IQGAP2
Other Names	Iqgap2;
Accession No.	Uniprot:Q13576
Calculated MW	181kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Background

Binds to activated CDC42 and RAC1 but does not seem to stimulate their GTPase activity. Associates with calmodulin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.