

TAF3 Rabbit mAb

Catalog No: #60198

Package Size: #60198-1 50ul #60198-2 100ul


Description

Product Name	TAF3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	TAF3 Antibody detects endogenous levels of total TAF3
Immunogen Description	A synthesized peptide derived from human TAF3
Other Names	TAF140; TAF3; TAFII140; TAFII155 protein;
Accession No.	Uniprot:Q5VWG9
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of TAF3 expression in HeLa cell lysate.

Background

Transcription factor TFIID is one of the general factors required for accurate and regulated initiation by RNA polymerase II. TFIID is a multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors. Required in complex with TBPL2 for the differentiation of myoblasts into myocytes. The complex replaces TFIID at specific promoters at an early stage in the differentiation process.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.