

UFC1 Rabbit mAb

Catalog No: #60209

Package Size: #60209-1 50ul #60209-2 100ul

Description

Product Name	UFC1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human
Specificity	UFC1 Antibody detects endogenous levels of total UFC1
Immunogen Description	A synthesized peptide derived from human UFC1
Other Names	HSPC155; Ubiquitin fold modifier conjugating enzyme 1; UFC1; Ufm1 conjugating enzyme 1;
Accession No.	Uniprot:Q9Y3C8
Calculated MW	19kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Background

E2-like enzyme which forms an intermediate with UFM1 via a thioester linkage.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.