

SEC24D Rabbit mAb

Catalog No: #60216

Package Size: #60216-1 50ul #60216-2 100ul

Description

Product Name	SEC24D Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	SEC24D Antibody detects endogenous levels of total SEC24D
Immunogen Description	A synthesized peptide derived from human SEC24D
Other Names	SEC24 family member D; SEC24D;
Accession No.	Uniprot:O94855
Calculated MW	113kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Background

Component of the COPII coat, that covers ER-derived vesicles involved in transport from the endoplasmic reticulum to the Golgi apparatus. COPII acts in the cytoplasm to promote the transport of secretory, plasma membrane, and vacuolar proteins from the endoplasmic reticulum to the Golgi complex.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.