

PWP1 Rabbit mAb

Catalog No: #60217

Package Size: #60217-1 50ul #60217-2 100ul

Description

Product Name	PWP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	PWP1 Antibody detects endogenous levels of total PWP1
Immunogen Description	A synthesized peptide derived from human PWP1
Other Names	PWP1;
Accession No.	Uniprot:Q13610
Calculated MW	75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Background

May play an important role in cell growth and/or transcription.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.