

NUP155 Rabbit mAb

Catalog No: #60227

Package Size: #60227-1 50ul #60227-2 100ul

Description

Product Name	NUP155 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	NUP155 Antibody detects endogenous levels of total NUP155
Immunogen Description	A synthesized peptide derived from human NUP155
Other Names	ATFB15; N155; Nup155;
Accession No.	Uniprot:O75694
Calculated MW	155kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Background

Essential component of nuclear pore complex. Nucleoporins may be involved both in binding and translocating proteins during nucleocytoplasmic transport.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.