

GNAQ Rabbit mAb

Catalog No: #60244

Package Size: #60244-1 50ul #60244-2 100ul

Description

Product Name	GNAQ Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	GNAQ Antibody detects endogenous levels of total GNAQ
Immunogen Description	A synthesized peptide derived from human GNAQ
Other Names	CMC1; G alpha Q; G protein alpha Q ; GAQ; GNAQ;
Accession No.	Uniprot:P50148
Calculated MW	42kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Background

Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.