

CSDE1 Rabbit mAb

Catalog No: #60246

Package Size: #60246-1 50ul #60246-2 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CSDE1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	CSDE1 Antibody detects endogenous levels of total CSDE1
Immunogen Description	A synthesized peptide derived from human CSDE1
Other Names	Csde1; NRU; UNR;
Accession No.	Uniprot:O75534
Calculated MW	89kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000 IP 1:50

Images



Western blot analysis of CSDE1 expression in (1) K562 cell lysate; (2) PC-12 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Background

RNA-binding protein. Required for internal initiation of translation of human rhinovirus RNA. May be involved in translationally coupled mRNA turnover.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.