

CRABP2 Rabbit mAb

Catalog No: #60251

Package Size: #60251-1 50ul #60251-2 100ul


Description

Product Name	CRABP2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	CRABP2 Antibody detects endogenous levels of total CRABP2
Immunogen Description	A synthesized peptide derived from human CRABP2
Other Names	Cellular retinoic acid binding protein 2; Cellular retinoic acid binding protein II ; Crabp2; CRABPII; RBP6;
Accession No.	Uniprot:P29373
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of CRABP2 expression in (1) MCF7 cell lysate; (2) Mouse skin lysate.

Background

Transports retinoic acid to the nucleus. Regulates the access of retinoic acid to the nuclear retinoic acid receptors.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.